


a major advance in cell counting and analyzing **introducing the new MODEL F Coulter Counter®**

Now, even easier, faster for the
most reliable and accurate counts

Featuring
New
Continuous
Twin Monitor*

1. Oscilloscope display
for complete
system check-out.
2. Orifice image
display for optical
aperture check-out.


*Patent applied for.



A completely
self-contained
compact unit
incorporating
new Coulter-
designed
transistorized
circuitry.

the second generation . . . MODEL F COULTER COUNTER

a direct development from the
world-famous Model A
Coulter Counter



Patented throughout the world.



Write or phone today for complete information or a demonstration in your own laboratory.

COULTER ELECTRONICS, INC.

590 West 20th Street, Hialeah, Florida Phone: 305-887-8131

BACTO - LATEX

0.81 MICRON

Product of *Dow Research & Difco Standardization*

¶ Bacto-Latex 0.81 micron is characterized by uniform particle size, batch reproducibility and biological inertness.

¶ An inert carrier for use in clinical and investigational tests including

RHEUMATOID ARTHRITIS

INFLAMMATORY DISEASES

TRICHINOSIS

LEPTOSPIROSIS

¶ Recommended for Rheumatoid Arthritis Tests of Singer and Plotz and modifications.

Literature available on request

DIFCO LABORATORIES DETROIT 1 MICHIGAN USA

BIOLOGICS

CULTURE MEDIA

REAGENTS

The logo consists of the word "DIFCO" in a bold, sans-serif font, enclosed within a stylized, dark, rounded rectangular border that resembles a film strip or a thick, dark frame.